

Hyperspectral Imaging: Bead Identification and Future Directions?

Emma Stavropoulos

Cissé Lab

Advisor: Dr. Jan-Hendrik Spille

Presentation 8/28/2017

Overall Goal: Simultaneously identify multiple fluorescent molecules across visible spectrum




Image of fluorescent beads





Image with spectral data

Uses a dispersive element (prism) in a secondary imaging pathway to disperse light


Requires a direct image and a dispersed (spectral) image to identify both position and color

□ Acquires all spectral information simultaneously, no extra filters or switching lasers!


Optical System



Optical System




Optical System




Nanogrid – 200 nm point sources

Direct




512 x 512 pixels


Spectral



1024 x 1024 pixels



Direct




Spectral

Note:

All images taken with filter Cube
with Semrock Dichroic Mirror
Reflects 405, 488, 561, 635 nm lasers


Nanogrid – 200 nm point sources

Direct




512 x 512 pixels

Spectral – 590 nm reference filter



1024 x 1024 pixels



Direct



Spectral


Direct

160 nm/pixel

Spectral


x ~ 30-40 nm/px

Deconvolve to disentangle spatial and spectral information




I. Transformation

Direct Path



Spectral Path:
590/10 nm reference




Spectral camera restriction


1. localize corresponding points on direct camera and spectral camera (590 nm)

I. Transformation

Direct Path




Spectral Path:
590/10 nm reference





Spectral camera restriction

1. localize corresponding points
on direct camera and spectral
camera (590 nm)


2. Affine transformation



II. Dispersion Map




Spectral image using 590 nm filter,
false colored orange



Composite image of 6 filters (465-710 nm),
false colored to match approximately

II. Dispersion Map


Spectral Path:
various filters





Spectral Path:
590/10 nm reference



As desired
linear dispersion ~ 3.6 nm/px!





Spectral image using 590 nm filter,
false colored orange



Composite image of 6 filters (465-710 nm),
false colored to match approximately


Mapping linear dispersion by imaging 7 filters and nanogrid




Three microbeads with single dyes

Goal:

Characterize three
beads on mixed slides
based on spectra




Direct Camera




Yellow-green bead (cyan) appears more shifted left than the red bead (yellow)


Direct camera – false colored composite




Spectral camera




Spectral camera intensity profiles



Direct camera intensity profiles



488 nm yellow-green bead spectrum




Imaging fixed tri-color beads

Direct Camera



Spectral Camera



642, 561, 488

642, 561, 488

642, 561, 488

642, 561, 488

Laser Lines (nm)

Dispersion:
Red ---> Blue spread






Imaging fixed tri-color beads

Direct Camera



Spectral Camera



642, 561, 488

642, 561, 488

642, 561, 488

642, 561, 488

Laser Lines (nm)


Dispersion:
Red ---> Blue spread



Future directions - Quantum Dots?

Pros

- Very narrow spectra
- Excitable by 405 nm laser
- Blink when isolated or can stably fluoresce
- Stable for a long time
- Bright!




Cons

- Much larger (tens of nm)
- Spectral blueing in oxygen environment and toxic to cells (mostly fixed)



Chan et al. Luminescent quantum dots for multiplexed biological detection and imaging. *Current opinion in biotechnology*. 2002; 13:40-46

Quantum dot pairs used as DNA ruler with superlocalization



0th order
(straight)

1st order
(deflected)



analytical
chemistry

Article
pubs.acs.org/ac


Superlocalization Spectral Imaging Microscopy of a Multicolor Quantum Dot Complex

Xingbo Shi,[†] Zhongqiu Xie,[†] Yuehong Song,[†] Yongjun Tan,[†] Edward S. Yeung,[†] and Hongwei Gai^{*,‡}


Quantum dot pairs used as DNA ruler with superlocalization

0th order
(straight)


1st order
(deflected)



1st order intensities



Localization model



analytical
chemistry

Article
pubs.acs.org/ac



Superlocalization Spectral Imaging Microscopy of a Multicolor Quantum Dot Complex

Xingbo Shi,[†] Zhongqiu Xie,[†] Yuehong Song,[†] Yongjun Tan,[†] Edward S. Yeung,[†] and Hongwei Gai^{*,‡}


Quantum dot pairs used as DNA ruler with superlocalization

0th order
(straight)

1st order
(deflected)



Localization model



analytical
chemistry

Article
pubs.acs.org/ac



Superlocalization Spectral Imaging Microscopy of a Multicolor Quantum Dot Complex

Xingbo Shi,[†] Zhongqiu Xie,[†] Yuehong Song,[†] Yongjun Tan,[†] Edward S. Yeung,[†] and Hongwei Gai^{*,‡}


Quantum dot pairs used as DNA ruler with superlocalization

0th order
(straight)

1st order
(deflected)



Localization model



analytical
chemistry

Article
pubs.acs.org/ac


Superlocalization Spectral Imaging Microscopy of a Multicolor Quantum Dot Complex

Xingbo Shi,[†] Zhongqiu Xie,[†] Yuehong Song,[†] Yongjun Tan,[†] Edward S. Yeung,[†] and Hongwei Gai^{*,‡}


Quantum dot pairs used as DNA ruler with superlocalization

0th order
(straight)


1st order
(deflected)



1st order intensities



Localization model



Conclusion:
estimated distance of 3 DNA lengths
within 2 nm of accepted value,
each with less than 7 nm uncertainty


analytical
chemistry

Article
pubs.acs.org/ac

Superlocalization Spectral Imaging Microscopy of a Multicolor Quantum Dot Complex


Xingbo Shi,[†] Zhongqiu Xie,[†] Yuehong Song,[†] Yongjun Tan,[†] Edward S. Yeung,[†] and Hongwei Gai^{*,‡}

Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)




460 + 650 QDs 1:1	
525 + 650 QDs 1:1	
590 + 650 QDs 1:1	
525 + 650 QDs 1:2	


Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)




Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)




Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)




Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)




Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)



Relative intensities determine ratio



Current hyperspectral setup may use similar QD pairs/triplets for **spectral barcoding** (...no superlocalization)



Hyperspectral imaging – what it (maybe) can and can't do

Requirements:

- fixed reference positions
for direct camera
- sparse but bright
- >50 ms exposure
- no blinking,
except for spectral barcodes
- no sub-pixel localization
experiments 😞

Hyperspectral imaging – what it (maybe) can and can't do


Requirements:

- fixed reference positions for direct camera
- sparse but bright
- >50 ms exposure
- no blinking, except for spectral barcodes
- no sub-pixel localization experiments 😞

Dynamics

Fixed objects with reference color (tagged in blue) and binding fluorescent probes (green/red)

- ☐ relative intensities reveal binding dynamics



Hyperspectral imaging – what it (maybe) can and can't do


Requirements:

- fixed reference positions for direct camera
- sparse but bright
- >50 ms exposure
- no blinking, except for spectral barcodes
- no sub-pixel localization experiments 😞

Dynamics

Fixed objects with reference color (tagged in blue) and binding fluorescent probes (green/red)


- relative intensities reveal binding dynamics



Multiplex Imaging

Blinking QD complexes bind to receptors on membrane

- frequency/type of cell receptors reveal information about cell



Thank you!

<3